

Beam Profiler



Features:

Spectral range: 300-1100nm

Beam diameter: 64 μm~22.4mm

Maximum power to 1W

PC software

SPECIFICATIONS

Product Name	CN0622VIS- B
Spectral Range (nm)	300-1100
Sensor Size (mm)	22.4*22.4
Pixel Size (μm)	6.4
Exposure Time (ms)	0.05~500
Shutter Type	Global Shutter cMOS
Beam Diameter	64μm~22.4mm
Maximum Power Density (W/cm ²)	10 (OD3.5)
Maximum Power (W)	1
Attenuator	OD0.3~OD5 (optional)
Communication Interface	10 GbE
External Triggering	Available
Dimensions	Reference structure diagram
Operating Environment	Temperature <0~35°C; Relative humidity <60%; Pointing test temperature: 15-25°C
Operating System	Windows 10/11, Optimum resolution: 1920*1080
Weight (g)	<500g
Warranty	Warranty

CN0622VIS- B



