



Beam Profiler



Features:

Spectral range:240-355nm

Beam diameter:1 μ m~100 μ m

Maximum power to 1W

PC software

SPECIFICATIONS:

Product Name	CN0205UV-XDT40
Spectral Range	240-355nm
Min focus	1 μ m
Max focus	100 μ m
Magnification	40x
Detection precision	$\pm 0.25 \mu$ m
Exposure Time	0.05ms~500ms
Maximum Power	1W
Attenuator	general band: OD1/OD2/OD3/OD4 (240-355nm)
Communication Interface	GigE, 支持 0.1POE
External Triggering	Available
Weight	<1.5Kg
Incoming light direction	Vertical (horizontal optional)
Defocusing range	-1mm, +1.3mm
Warranty	1 year

CN0205UV-XDT40

