



## Beam Profiler



### Features:

Spectral range:400-800nm

Beam diameter:5 $\mu$ m~500 $\mu$ m

Maximum power to 1W

PC software

### SPECIFICATIONS:

Product Name	CN0205VIS-XV
Spectral Range	400-800nm
Min focus	5 $\mu$ m
Max focus	500 $\mu$ m
Magnification	10x
Detection precision	$\pm 1\mu$ m
Exposure Time	0.05ms~500ms
Maximum Power	1W
Attenuator	general band: OD1/OD2/OD3/OD4 (400-800nm)
Communication Interface	USB3.0
External Triggering	Available
Weight	<1.5Kg
Incoming light direction	Vertical (horizontal optional)
Defocusing range	-1mm, +1.3mm
Warranty	1 year

## CN0205VIS-XV

