



Beam Profiler



Features:

Spectral range:900-1100nm

Beam diameter:2.5 μ m~250um

Maximum power to 3W

PC software

SPECIFICATIONS:

Product Name	CN0307VIS-XN20
Spectral Range	900-1100nm
Min focus	2.5 μ m
Max focus	250um
Magnification	20x
Detection precision	$\pm 1 \mu m$
Exposure Time	0.05ms~500ms
Maximum Power	3W
Attenuator	general band: OD1/OD2/OD3/OD4
Communication Interface	USB3.0
External Triggering	Available
Incoming light direction	Vertical (horizontal optional)
Defocusing range	-1mm, +1.3mm
Weight	<1.5Kg
Warranty	1 year

CN0307VIS-XN20

